

# INFANT WARMER CARE UNIT

## KEY FEATURES

- Fully Automatic Heater Controller with Baby Skin Sensor
- Adjustable Heat Setting up to 37 C
- Heater Active/Deactivate push (on/off)
- Automatic Heater Power (0 to 100%) Control Heating System
- On Screen Heater Power / Real Temp / Set Temp Readings
- High / Low Heating / Sensor Break Alarms
- Skin Sensor accidental Drop Down Alarm
- YSI400 Standard Skin Sensor with High accuracy (+- 0.25C @ 25.0C)
- Faults Indication Onscreen and Sound Buzzer
- Alarm Buzzer Mute
- Examination Lamp / Left / Right Tubes
- Heater Power Solid State Control no Mechanical Relay or Switch
- IV Stand / Cardio Monitor Shelf
- High Grade Powder Coated MS Sheets Design
- One Drawer and One Cabinet at lower bed
- Bed Up / Down Manual Adjustment (-5 to +10 degree)
- Touch Push and Splash Prof Front Panel



NOTE:  
Any Modification Available on Demand

**MODEL: BBL-7SHTD**



**ALLIED TECHNO CONTROL SYSTEMS**  
JHUMRA ROAD, MANSOORABAD,  
FAISALABAD, PAKISTAN-38000.  
+92-321-6622966  
ALLIEDTECHNOLOGIST@GMAIL.COM

www.alliedtechno.com

## SPECIFICATIONS

<b>PLC Controller Series</b>	7S HTD Series Warmer Controller
<b>Dimensions</b>	72 H x 26 W x 32 D (inch)
<b>Net Weight</b>	50Kg (approx.)
<b>Connection Type</b>	PAK Standard UL Certified
<b>Power</b>	200 to 260 VAC - 1050 Watts max. 50-60Hz
<b>Heater Type</b>	INFRA RED HEATING SYSTEM
<b>Operating Environment</b>	5 to 40°C , humidity: 10 to 80% RH, non-condensing
<b>Enclosure / Body / Frame</b>	Flame-Retardant Metal Sheet - ABS material(UL-94) - Acrylic
<b>Warm-up Time</b>	Depend On Environment Temp - 1C Per 2 Minutes Maximum @ 25C Ambient
<b>Display Resolution</b>	0.1°C (with YSI400 Skin Probe RTD)
<b>Wheel Type</b>	Rubber Wheel 4 inch ( 2x Free / 2x Brake)
<b>Side Bumper Protection</b>	N / A
<b>Paint Type</b>	High Grade Epoxy Powder Coating (Customized Available)
<b>Side Lights</b>	Standard White Light ( Therapy Light Optional)
<b>Examination Light</b>	Natural Warm Light (Skin Hue Compatible)